

StellarNet Spectrometer Python Module

This module requires Python 2.7 and depends on pyusb 1.0.x. One way to install it is using pip:

Install pip: <http://www.pip-installer.org/en/latest/installing.html>
Install pyusb: `pip install pyusb==1.0.0a3`

After plugging in a device, you should be able to test the basic operation of the driver:

```
sudo python stellarnet.py info
```

Access to USB devices on Linux systems is restricted and requires root access. However, it is possible to overcome this restriction by installing the custom udev rules file provided: 99-local.rules. Install this file as /etc/udev/rules.d/99-local.rules and reboot the system to provide non-root access to StellarNet devices.

Class description

class stellarnet.**StellarNet**(*device*)¶

Represents a StellarNet spectrometer.

COEFF_C1_ADDR = 128¶

The address of the string containing coefficient C1

COEFF_C2_ADDR = 160¶

The address of the string containing coefficient C2

COEFF_C3_ADDR = 192¶

The address of the string containing coefficient C3

COEFF_C4_ADDR = 224¶

The address of the string containing coefficient C4

DEVICE_ID_ADDR = 32¶

The address of the stored string containing device identification.

Exceptions

exception stellarnet.**ArgRangeError**(*message*)¶

Raised when argument is out of range.

exception stellarnet.**ArgTypeError**(*message*)¶

Raised when argument type is incorrect.

exception stellarnet.**ArgumentError**(*message*)¶

Raised when argument in error.

exception stellarnet.**NotFoundError**(*message*)¶

Raised when USB device cannot be found.

A typical application will involve an operational sequence similar to the following:

```
find_devices()
```

```
get_device_id()
```

```
get_config()
```

```
set_config(**kwargs)
```

```
read_spectrum()
```

Alphabetic Listing

compute_lambda(*pixel*)¶

Compute lambda from the pixel index. Returns the pixel's wavelength (float).

Parameters:

pixel – Integer; the pixel index on which to perform the computation.

get_config()¶

Gets the device configuration. Returns a dictionary.

get_device_id()¶

Gets the device id. Returns a string.

get_stored_bytes(*address*)¶

Get stored bytes. Returns bytearray.

Parameters:

address – Integer; the address of the string to get.

get_stored_string(*address*)¶

Get stored bytes. Returns string.

Parameters:

address – Integer; the address of the string to get.

print_info()¶

Print device information.

read_spectrum()¶

Reads and returns a spectrum from the spectrometer. Returns an array of short integers.

See [stellarNet.set_config\(\)](#) for a description of the parameters that control the operation of the spectrometer or the post-processing of the spectrum.

set_config(*kwargs*)**¶

Sets the device configuration.

Parameters:

- **int_time** – (optional) Integer; the integration time in milliseconds.
- **x_timing** – (optional) Integer; the XTiming rate.
- **x_smooth** – (optional) Integer; the boxcar smoothing window size.
- **scans_to_avg** – (optional) Integer; the # of scans to be averaged together.

- **temp_comp** – (optional) Integer; temperature compensation (not implemented).

set_stored_bytes(*address*, *data*)¶

Set stored bytes.

Parameters:

- **address** – Integer; the address of the string to set.
- **data** – String; the string value to be set

exception `stellarnet.StellarNetError`(*message*)¶

Base class for StellarNet errors.

exception `stellarnet.TimeoutError`(*message*)¶

Raised when device operation times out.

find_devices()¶

Find all USB-connected StellarNet devices.

This function returns a tuple of StellarNet objects or raises `NotFoundError` if no devices are found.

main(*argv=None*)¶

Exercise StellarNet spectrometer driver.